



		OP-S6	OP-S7	OP-S8	OP-S10	OP-S14
Up/down movement		Electrical (motor)	Electrical (motor)	Electrical (motor)	Electrical (motor)	Electrical (motor)
Reclining movement		Manual (gas)	Electrical (motor)	Electrical (motor)	Electrical (motor)	Electrical (motor)
Movement control	Foot Pedal	Yes	Yes	Yes	Yes (Optional)	Yes (Optional)
	Keyboard	--	--	2 (Optional)	4	4
Memories		--	2	2	6	6
Rotation Brake		Manual	Manual	Manual	Electro-magnetic	Electro-magnetic
Fold down armrest (hide away)		Yes	Yes	Yes	Yes	Yes
Rotation		360°	370°	300°	300°	300°
Trendelenburg Position		Yes	Yes	Yes	Yes	Yes
Synchronized sliding backrest		--	--	--	Yes	Yes
80° Backrest position		--	--	--	Yes	Yes
TFT Monitor		--	--	--	--	Optional
Monitor with integrated camera		--	--	--	--	Optional
LED Light source		--	--	--	--	Optional
Examination Spotlight		--	--	--	--	Optional

Examination Chairs